


Why Do You Need Preortho

- ▶ Malocclusions
- ▶ Remove by oneself
- ▶ Preform type
- ▶ Utilized for MFT

Preortho improves habits though MFT.

Preortho provides you to proper alignment though MFT.

There are maximum results prior to permanent dentition.



» Contact Us

Authorized Distributor



+82-51-556-2070
+82-10-4032-8886
optimaortho.korea@gmail.com
www.ipreortho.com




YouTube
Preortho Official Channel
<https://www.youtube.com/@Preorthoofficial>



Preortho

Pre-Orthodontic Treatment & MyoFunctional Therapy Device



www.ipreortho.com




Preortho 3 Features

1. Soft & Comfortable

- Comfortability is Excellent
- Durable Appliance

1 Year
Warranty



2. Pre-Formed Appliance

- Hassle Free Appliance
- Minimum Adjustment for Customized Fit
- Mouth Piece Type of Orthodontic Appliance

3. For MFT

- Better Nasal Breathing vs Mouth Breathing
- Improves Speech Impediment & Swallowing



Inventor

Dr. Otsuka Atsushi


- Tokyo Dental College
- DDS, MD, Orthodontist
- Otsuka Orthodontic Clinic
- Inventor of Preortho
- Preortho Chief Clinical Consultant




Preortho Product List

	Initial	S size	M size	L size
Type I Maxillary protrusion Crowding Deep bite Retention	Initial	S	M	L
		 Hard Type Soft Type	 Hard Type Soft Type	 Hard Type
Type II Open bite Retention	—	—	M  Hard Type	—
Type III Cross bite Retention	Initial COMING SOON	S  Hard Type Soft Type	M  Hard Type	L  Hard Type

>> Type I



>> Type II



>> Reshape of Preortho

Preortho warms the device to expand to fit the patient, allowing for partial fabrication (customize) with a burr or surgical scissors.



Preortho Official Channel
"How to expand Preortho Appliance"
<https://www.youtube.com/watch?v=Ar6iVt-JGHo>

>> Type III

Initial



COMING
SOON

S



60.8mm

M



62.4mm


L



68.3mm


>> Buccinator mechanism

Rear view



※The illustration is Type I.

Side view



※The illustration is Type I.

Provide a normal overjet.

Difficult to achieve incisal contact during the treatment process.

>> A Message From The Inventor

As an orthodontic specialist with nearly 30 years of clinical experience, I have treated numerous cases of pediatric malocclusion using devices such as FKO, Bimler, and recently developed mouthpiece-type appliances and multi-bracket appliances.

About a decade ago, I began working on developing a treatment method that is less burdensome and safer for patients. After many research and refinement, the current Preortho system for pediatric orthodontics was established. In recent years, we have successfully developed both the treatment protocol and the appliance, making it suitable for routine clinical application. Many dental specialists have adopted this device and treatment method in their practices, resulting in numerous successful outcomes.

Moving forward, I am dedicated to further contributing to the field by sharing valuable clinical experiences and academic resources.



Inventor
Dr. Atsushi Otsuka

- Tokyo Dental College
- DDS, MD, Orthodontist
- Otsuka Orthodontic Clinic
- Inventor of Preortho
- Preortho Chief Clinical Consultant

>> Type II

Before



▶ ▶ ▶ ▶ ▶ ▶ ▶ ▶ ▶ ▶

23months later

Age: 7yrs old, Male

Open bite



Preortho Only

After



by Ulsan BS Orthodontics Clinic / Dr. YoungGyu Lee

>> Type III

Before



▶ ▶ ▶ ▶ ▶ ▶ ▶ ▶ ▶ ▶

22months later

Age: 8yrs old, Male

Anterior Crossbite



Preortho Only

After



by Ulsan BS Orthodontics Clinic / Dr. YoungGyu Lee

>> Type I + Type II

Before



▶ ▶ ▶ ▶ ▶ ▶ ▶ ▶ ▶ ▶

16months later

Age: 5yrs old, Male

Overbite protrusion

Preortho Only T1-SS

After



Before



▶ ▶ ▶ ▶ ▶ ▶ ▶ ▶ ▶ ▶

12months later

Age: 6yrs old, Male

Open bite

Preortho Only T2-MH

After



by Yonsei Dreaming Dental Clinic / Dr. HyunMook Lim

Before



▶ ▶ ▶ ▶ ▶ ▶ ▶ ▶ ▶ ▶

40months later

Age: 6yrs old, Male

False Crossbite



Preortho Only

After



by Yonsei Dreaming Dental Clinic / Dr. KiHoon Joo



The Preortho treatment Flow Chart

The supplier will provide the Preortho treatment Flow Chart for your reference.

If you have any clinical question, please find the user group on Facebook.

Preortho User Group

QR Code ▶▶▶

<https://m.site.naver.com/0wS5d>

